


低 VOC & 低アウトガス

*LaserClean*TM

光学ブレットボード



- 真空対応製品のクリーン度を低価格で提供
- 独自の製造 & クリーニング処理で低アウトガスを保証
- サイエнтиフィックグレード性能
- シロキサン等の VOC、有害物質を低減:
- アルマイト加工 SA2 アルミ製もラインナップ
- レーザー & メディカル分野向け

製品一覧

モデル	概要
LC-25-2	Optical Breadboard, LaserClean, 4.8 mm Skin, 24 x 60 x 2.3 in., 1/4-20
LC-25-4	Optical Breadboard, LaserClean, 4.8 mm Skin, 24 x 60 x 4.3 in., 1/4-20
LC-35-2	Optical Breadboard, LaserClean, 4.8 mm Skin, 36 x 60 x 2.3 in., 1/4-20
LC-35-4	Optical Breadboard, LaserClean, 4.8 mm Skin, 36 x 60 x 4.3 in., 1/4-20
M-LC-25-2	Optical Breadboard, LaserClean, 4.8 mm Skin, 600 x 1500 x 59 mm, M6
M-LC-25-4	Optical Breadboard, LaserClean, 4.8 mm Skin, 600 x 1500 x 110 mm, M6
M-LC-35-2	Optical Breadboard, LaserClean, 4.8 mm Skin, 900 x 1500 x 59 mm, M6
M-LC-35-4	Optical Breadboard, LaserClean, 4.8 mm Skin, 900 x 1500 x 110 mm, M6
M-SA2-06x06-LC	Solid Aluminum Optical Breadboard, LaserClean, 150 x 150 mm, M6
M-SA2-08x08-LC	Solid Aluminum Optical Breadboard, LaserClean, 200 x 200 mm, M6
M-SA2-08x10-LC	Solid Aluminum Optical Breadboard, LaserClean, 250 x 200 mm, M6
M-SA2-11-LC	Solid Aluminum Optical Breadboard, LaserClean, 300 x 300 mm, M6
M-SA2-12-LC	Solid Aluminum Optical Breadboard, LaserClean, 600 x 300 mm, M6
M-SA2-12x18-LC	Solid Aluminum Optical Breadboard, LaserClean, 450 x 300 mm, M6
M-SA2-13-LC	Solid Aluminum Optical Breadboard, LaserClean, 900 x 300 mm, M6
M-SA2-18x18-LC	Solid Aluminum Optical Breadboard, LaserClean, 450 x 450 mm, M6
M-SA2-18x24-LC	Solid Aluminum Optical Breadboard, LaserClean, 600 x 450 mm, M6
M-SA2-18x30-LC	Solid Aluminum Optical Breadboard, LaserClean, 750 x 450 mm, M6
M-SA2-22-LC	Solid Aluminum Optical Breadboard, LaserClean, 600 x 600 mm, M6
M-SA2-23-LC	Solid Aluminum Optical Breadboard, LaserClean, 900 x 600 mm, M6
M-SA2-30x30-LC	Solid Aluminum Optical Breadboard, LaserClean, 750 x 750 mm, M6



株式会社日本レーザー

東京本社 〒169-0051 東京都新宿区西早稲田 2-14-1 TEL 03-5285-0853
 大阪支店 〒533-0033 大阪市東淀川区東中島 1-20-12 TEL 06-6323-7286
 名古屋支店 〒460-0003 名古屋市中区錦 3-1-30 TEL 052-205-9711
 e-mail: newport@japanlaser.jp / http://www.newport-japan.jp