

狭い間隔で光学部品を取付け。コンパクトな光学系の構築に。

## 2倍穴密度 アルミブレッドボード・プレート



- 取付け穴配列の密度が標準の2倍
- M6 or 1/4-20 取付け穴
- 12.7mm 厚アルミ製
- 作業面: 低反射アルマイト処理

### 製品一覧

モデル	概要
M-SA2-06x06-D	Solid Aluminum Optical Breadboard, 150 x 150 mm, Double Density M6 Grid
M-SA2-06x12-D	Solid Aluminum Optical Breadboard, 150 x 300 mm, Double Density M6 Grid
M-SA2-06x18-D	Solid Aluminum Optical Breadboard, 150 x 450 mm, Double Density M6 Grid
M-SA2-06x24-D	Solid Aluminum Optical Breadboard, 150 x 600 mm, Double Density M6 Grid
M-SA2-08x08-D	Solid Aluminum Optical Breadboard, 200 x 200 mm, Double Density M6 Grid
M-SA2-11-D	Solid Aluminum Optical Breadboard, 300 x 300 mm, Double Density M6 Grid
M-SA2-12-D	Solid Aluminum Optical Breadboard, 300 x 600 mm, Double Density M6 Grid
M-SA2-12x18-D	Solid Aluminum Optical Breadboard, 300 x 450 mm, Double Density M6 Grid
M-SA2-18x18-D	Solid Aluminum Optical Breadboard, 450 x 450 mm, Double Density M6 Grid
M-SA2-22-D	Solid Aluminum Optical Breadboard, 600 x 600 mm, Double Density M6 Grid
SA2-06x06-D	Solid Aluminum Optical Breadboard, 6 x 6 in., Double Density 1/4-20 Grid
SA2-06x12-D	Solid Aluminum Optical Breadboard, 6 x 12 in., Double Density 1/4-20 Grid
SA2-06x18-D	Solid Aluminum Optical Breadboard, 6 x 18 in., Double Density 1/4-20 Grid
SA2-06x24-D	Solid Aluminum Optical Breadboard, 6 x 24 in., Double Density 1/4-20 Grid
SA2-08x08-D	Solid Aluminum Optical Breadboard, 8 x 8 in., Double Density 1/4-20 Grid
SA2-11-D	Solid Aluminum Optical Breadboard, 12 x 12 in., Double Density 1/4-20 Grid
SA2-12-D	Solid Aluminum Optical Breadboard, 12 x 24 in., Double Density 1/4-20 Grid
SA2-12x18-D	Solid Aluminum Optical Breadboard, 12 x 18 in., Double Density 1/4-20 Grid
SA2-18x18-D	Solid Aluminum Optical Breadboard, 18 x 18 in., Double Density 1/4-20 Grid
SA2-22-D	Solid Aluminum Optical Breadboard, 24 x 24 in., Double Density 1/4-20 Grid



株式会社日本レーザー

東京本社 〒169-0051 東京都新宿区西早稲田 2-14-1 TEL 03-5285-0853  
 大阪支店 〒533-0033 大阪市東淀川区東中島 1-20-12 TEL 06-6323-7286  
 名古屋支店 〒460-0003 名古屋市中区錦 3-1-30 TEL 052-205-9711  
 e-mail: [newport@japanlaser.jp](mailto:newport@japanlaser.jp) / <http://www.newport-japan.jp>